EN	ENGLISH
----	---------

Trypan Blue Ophthalmic Solution 0.06%

INTRODUCTION

S_{0x}-BLUE is a sterile solution of trypan blue (an acid di-azo group dye), a selective tissue-staining agent used as a medical aid in ophthalmic surgery.

Molecular formula: C₂₄H₂₄N₄Na₂O₄S₂
Molecular weight: 960.8
Chemical name: 3,3'- [(3,3'-dimethyl-4,4'-biphenylene) bis (azo)] bis(5-amino-4-hydroxy-2,7-naphthalenedisulfonic acid) tetra sodium salt
Composition per 1ml: Trypan Blue 0.6mg; sterile aqueous base q.s.

S_{0x}-BLUE is available either in 1ml vial (REF. 84000) or in 1ml glass pre-filled syringe accompanied by a 23G cannula (REF. 84001).

INTENDED USE

In cataract surgery, it is used to stain the anterior capsule of the human lens. As an inadequate visualization of the anterior capsule may result in improper or incomplete capsulorhexis carrying a high risk of radial capsular tears and other associated complications, the staining of the capsule improves the visualization and facilitates the performance of Capsulorhexis. More specifically, during cataract surgery:

- ✓ It offers enhanced visualization of the anterior capsule in eyes with mature cataract or narrow pupil
- ✓ It provides a persistent and clear outline of the peripheral rim of the capsulorhexis by contrast between stained rim and adjacent lens mass
- ✓ It reduces the risk of capsulorhexis-related complications by better visualization of radial capsule tears
- ✓ It provides visualization of the "lost" edge of the capsulorhexis.

CONTRAINDICATIONS

S_{0x}-BLUE should not be used when a non-hydrated (dry state) hydrophilic acrylic intraocular lenses (IOL) is planned to be implanted because the dye may be absorbed by the IOL and stain it.

WARNINGS / PRECAUTIONS

- Read instructions before use
- Inspect and check the product before use. Don't use it, if the packing is damaged
- The cannula and the syringe are for single use only
- For intraocular use only, not for injection.
- Do not use the product after its expiry date
- Not to be used on children or pregnant women
- Product should not be re-sterilized
- Mishandling or improper transport can cause damage to the device or packaging
- The product is non-toxic, sterile and non-pyrogenic, if packaging is intact
- Dosage as directed by the physician.

ADVERSE REACTIONS

Discoloration of specific intraocular lenses (see *Contraindications*) may occur, as well as an inadvertent staining of the posterior lens capsule and vitreous face, which is generally self-limited, lasting up to one week.

INSTRUCTIONS FOR USE

- Store in a cool and dry place (10-30 °C)
- Check package integrity before use, to ensure that the product has remained sterile
- Contains no preservative. If solution becomes dark brown or floating particles found, it should be discarded
- Each unit of S_{0x}-BLUE is sterilized, hence does not need autoclaving before surgery
- The normal precautions associated with ocular microsurgery should be observed
- The quantity injected into the anterior chamber of the eye must be adjusted according to the volume of the aqueous humour and the anatomical structure to be protected
- Remove the entire product by irrigation and/or aspiration at the end of the procedure; mechanical blockage of drainage at the trabecular level may cause a transient increase in intra-ocular pressure after surgery
- All post-operative inflammatory reactions (iritis, hypopyon, uveitis) are inherent in surgical procedures involving the anterior chamber of the eye, and no relationship with the product has been established
- Must not be reused. Otherwise, it may cause serious problems to patients.

The manufacturer does not hold any responsibility for any consequences deriving from improper use of the product.

DE	DEUTSCH
----	---------

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